

Short Scan Using XCIST

Botao Zhao

October 10, 2025

1 Introduction

This document explains how to perform short-scan reconstruction using XCIST. XCIST supports two reconstruction modes: `fdk_equiAngle` and `helical_equiAngle`. `fdk_equiAngle` is used for axial scans, where the patient table remains stationary (table speed = 0). `helical_equiAngle` is used for helical scans, where the patient table moves continuously (table speed $\neq 0$). The table speed can be configured in the `protocol.cfg` file, while the reconstruction type is specified in the `recon.cfg` file.

```
1 protocol.scanTypes = [1, 1, 1, 1] # flags for airscan, offset scan, phantom scan, prep
2
3 # Table and gantry
4 protocol.scanTrajectory = "Gantry_Helical" # name of the function that defines the scanning trajectory and model
5 protocol.viewsPerRotation = 500 # total numbers of view per rotation
6 protocol.viewCount = 500 # total number of views in scan
7 protocol.startViewId = 0 # index of the first view in the scan
8 protocol.stopViewId = protocol.startViewId + protocol.viewCount - 1 # index of the last view in the scan
9 protocol.airViewCount = 1 # number of views averaged for air scan
10 protocol.offsetViewCount = 1 # number of views averaged for offset scan
11 protocol.rotationTime = 1.0 # gantry rotation period (in seconds)
12 protocol.rotationDirection = 1 # gantry rotation direction (1=CCW, -1 C CW, seen from table foot-end)
13 protocol.startAngle = 0 # relative to vertical y-axis (n degrees)
14 protocol.tableSpeed = 0 # speed of table translation along positive z-axis (in mm/sec)
15 protocol.startZ = 0 # start z-position of table
16 protocol.tiltAngle = 0 # gantry tilt angle towards negative z-axis (in degrees)
17 protocol.wobbleDistance = 0.0 # focalspot wobble distance
18 protocol.focalspotOffset = [0, 0, 0] # focalspot position offset
19
20 # X-ray tube technique and filtration
21 protocol.mA = 500 # tube current (in mA)
22 protocol.spectrumCallback = "Spectrum" # name of function that reads and models the X-ray spectrum
23 protocol.spectrumFilename = "xcist_kv120_tar7_bin1.dat" # name of the spectrum file
24 protocol.spectrumUnit_mm = 1; # Is the spectrum file in units of photons/sec/mm^2/<current>?
25 protocol.spectrumUnit_mA = 1; # Is the spectrum file in units of photons/sec/<area>/mA?
26 protocol.spectrumScaling = 1; # scaling factor, works for both mono- and poly-chromatic spectra
27 protocol.bowtie = "Large.txt" # name of the bowtie file (or [] for no bowtie)
28 protocol.filterCallback = "Xray_Filter" # name of function to compute additional filtration
29 protocol.flatFilter = ['air', 0.001] # additional filtration - materials and thicknesses (in mm)
30 protocol.dutyRatio = 1.0 # tube ON time fraction (for pulsed tubes)
31 protocol.maxPrep = -1 # set the upper limit of prep, non-positive will disable this feature
```

```
1 recon.fov = 500.0 # diameter of the reconstruction field-of-view (in mm)
2 recon.imageSize = 512 # number of columns and rows to be reconstructed (square)
3 recon.sliceCount = 1 # number of slices to reconstruct
4
5 recon.sliceThickness = 0.579 # reconstruction slice thickness AND inter-slice interval (in mm)
6 recon.centerOffset = [0.0, 0.0, 0.0] # reconstruction offset relative to center of rotation (in mm)
7 recon.reconType = 'fdk_equiAngle' # Name of the recon function to call
8 recon.kernelType = 'standard' # 'R-L' for the Ramachandran-Lakshminarayanan (R-L) filter, rectangular window
9 # 'S-L' for the Shepp-Logan (S-L) filter, sinc window function
10 # 'soft', 'standard', 'bone' for kernels similar to those on clinical scanners
11 recon.startAngle = 0 # in degrees; 0 is with the X-ray source at the top
12 recon.unit = 'HU' # '/mm', '/cm', or 'HU'
13 recon.mu = 0.02 # in /mm; typically around 0.02/mm
14 recon.huOffset = -1000 # unit is HU, -1000 HU by definition but sometimes something else is present
15 recon.printReconParameters = False # Flag to print the recon parameters
16 recon.saveImageVolume = True # Flag to save recon results as one big file
17 recon.saveSingleImages = False # Flag to save recon results as individual images
18 recon.displayImagePictures = False # Flag to display the recon results as .png images
19 recon.saveImagePictureFiles = False # Flag to save the recon results as .png images
20 recon.displayImagePictureAxes = False # Flag to display the axes on the .png images
21 recon.displayImagePictureTitles = False # Flag to display the titles on the .png images
22
```

2 Running Short Scan

The easiest way to do short scan is setting the value of `protocol.viewCount` to be less than `protocol.viewsPerRotation`, in `protocol.cfg` file. For example, if you set `protocol.viewsPerRotation` = 100 and `protocol.viewCount` = 50, that's equivalent to "half scan".

```
1 protocol.scanTypes = [1, 1, 1, 1] # flags for aiscan, offset scan, phantom scan, prep
2
3 # Table and gantry
4 protocol.scanTrajectory = "Gantry_Helical" # name of the function that defines the scanning trajectory and model
5 protocol.viewsPerRotation = 100 # total numbers of view per rotation
6 protocol.viewCount = 50 # total number of views in scan
7 protocol.startViewId = 0 # index of the first view in the scan
8 protocol.stopViewId = protocol.startViewId + protocol.viewCount - 1 # index of the last view in the scan
```


2.1 Running Short Scan - Axial Scanning

To run short scan with axial scanning mode, you need to set `recon.reconType = 'fdk_equiAngle'` in `recon.cfg`. However, the default code doesn't allow you to do short scan in axial scan mode (It will gives you index error). So we have to modify the source code `fdk_equiAngle.py` in the XCIST package.

Run `pip show gecatsim` in Anaconda Prompt:

```
Anaconda Prompt
(base) C:\Users\Botao>pip show gecatsim
Name: gecatsim
Version: 1.5.4
Summary: Simulation toolkit for X-ray based cancer imaging
Home-page: https://github.com/xcist/documentation/wiki
Author: Mingye Wu, Paul FitzGerald, James Jobin, Jiayong Zhang, Nandan Reddy, Anupama Debnath, Brion Sarachan, Bruno De Man
Author-email: Mingye.Wu@gehealthcare.com
License: BSD 3-Clause License
Location: C:\Users\Botao\anaconda3\Lib\site-packages
Requires: matplotlib, numpy, scipy, tqdm
Required-by:
```

This should give you the location of where the package is installed at. Then go to `gecatsim\reconstruction\pyfiles`. Open `fdk_equiAngle.py`.



C:\Users\Botao\anaconda3\Lib\site-packages\gecatsim\reconstruction\pyfiles			
Name	Date modified	Type	Size
__pycache__	10/7/2025 9:43 AM	File folder	
createHSP.py	9/20/2025 4:09 PM	PY File	3 KB
fdk_equiAngle.py	9/20/2025 4:09 PM	PY File	10 KB
helical_equiAngle.py	9/20/2025 4:09 PM	PY File	11 KB
mapConfigVariablesToFDK.py	9/20/2025 4:09 PM	PY File	2 KB
mapConfigVariablesToHelical.py	9/20/2025 4:09 PM	PY File	3 KB
recon.py	9/20/2025 4:09 PM	PY File	7 KB

At line 138, replace `Dg = Dgy` with `Dg = np.append(Dgy, np.zeros((YL, ZL, ProjScale - cfg.protocol.viewCount)), axis=2)`:


```

129 ##### pre-weighting for ramp-filter
130
131
132 print(" Pre-weighting the filter...")
133 for Yindex in range(YL):
134     for zindex in range(ZL):
135         Dgy[Yindex, zindex, :] = (DistD / np.sqrt(DistD ** 2 + ((zindex - ZCtr) * DeltaZ) ** 2)) * ProjData[Yindex,
136             zindex, :] * math.cos((Yindex - YCtr) * DeltaUW)
137
138 #Dg=Dgy
139 Dg = np.append(Dgy, np.zeros((YL, ZL, ProjScale - cfg.protocol.viewCount)), axis=2) # try padding (cut unwanted section from the
140 scanned data)
141

```

What this does is it pads the missing section from a short scan (pretending it's a full scan but fill the missing section with 0), and then reconstructs using pre-defined view counts.

Running the main script with `protocol.viewsPerRotation = 100` and `protocol.viewCount = 50`. You should get the result shown below:



2.2 Running Short Scan - Helical Scanning

Running Short Scan with Helical Scanning mode is much easier. Simply set `recon.reconType = 'helical_equiAngle'` in `recon.cfg`. You don't need to modify anything, but the table speed must be non-zero otherwise the recon result will be pure dark image.